



VWR Storage & transport containers




Thermal protection for
Perishables & Valuables

Refrigerant Bottles

Refrigerant Bricks

Containers, Shipping,
Insulated PUR

Refrigerant Gel Packs



Refrigerant Bottles

- Blow-moulded rigid plastic bottles
- Durable gels are multi-purpose to withstand rugged packing environments
- Bottles are reusable

Rigid Bottle refrigerant gels are designed to provide reliable temperature protection for an array of products and packages. These gels deliver the stability you require from your refrigerant. The bottles are strong for rough transport and handling, reusable, and return to a gel state after thawing.

Description	Dimensions (hwxwd) (mm)	Qty/Pk	Cat. No.
Cold Bottle	228 x 101 x 51	36	216-0186
Cold Bottle	216 x 177 x 51	24	216-0187



Refrigerant Bricks

- Semi-rigid foam that never changes shape
- Maintains designated packing configuration throughout trip
- Vacuum-formed to ensure an exact release of thermal energy

Foam bricks are the ideal refrigerant for when space is limited. Bricks never change shape during freezing, thawing, or transit. These refrigerants are available in a variety of sizes. Each brick is vacuum sealed for consistent weight and shape.

Description	Dimensions (hwxwd) (mm)	Weight (g)	Qty/Pk	Cat. No.
Cold Brick	114 x 102 x 38	426	36	216-0188
Cold Brick	229 x 102 x 38	880	18	216-0189
Cold Brick	286 x 241 x 25	1590	12	216-0190
Cold Brick	114 x 102 x 38	199	72	216-0191



Containers, Shipping, Insulated PUR

- Polyurethane insulated containers provide extensive protection for long range journeys
- Durable, one-piece closed-cell polyurethane foam insulated shippers are CFC and HCFC free. They also comply with international clean air regulations
- Ensures frozen or refrigerated products from moderate and extreme ambient temperature exposure

These high performance insulated **Polyurethane** containers meet the requirements of the most demanding temperature sensitive contents. Polyurethane insulation offers greater insulating power. Thick walls offer protective shock absorption and cushioning during transportation. Containers include foam plug to seal the top, cushion, insulate and vent dry ice.

Description	Dimensions (hwxwd) (mm)	Qty/Pk	Cat. No.
Shipping Container	279 x 203 x 177	3	216-0211
Shipping Container	413 x 235 x 260	2	216-0212
Shipping Container	470 x 464 x 425	1	216-0213
Shipping Container	483 x 400 x 304	1	216-0214



Refrigerant Gel Packs

- Non-toxic protection for shipping temperature sensitive materials
- Phase changes at 0°C/+32°F
- Returns to gel state after thawing

As an alternative to wet ice, gel pack refrigerants provide proven quality and reliability to deliver temperature sensitive packages in top condition. Gel packs are reusable and available in a range of sizes.



Description	Dimensions (hxwxh) (mm)	Weight (g)	Qty/Pk	Cat. No.
Gel Pack	165 x 89 x 19	227	72	216-0192
Gel Pack	153 x 140 x 25	341	48	216-0193
Gel Pack	152 x 152 x 25	454	36	216-0194
Gel Pack	230 x 140 x 32	682	24	216-0195
Gel Pack	267 x 140 x 32	909	18	216-0196
Gel Pack	100 x 150 x 20	-	72	216-0215
Gel Pack	175 x 150 x 25	-	36	216-0216
Gel Pack	200 x 150 x 32	-	24	216-0217

Containers, Shipping, Insulated EPS

- Expanded polystyrene insulated shipping containers are ideal for the safe storage and delivery of refrigerated, frozen and other perishable products
- Choose EPS or EPS with corrugated carton
- Economical, lightweight, reusable, recyclable



Insulated shipping containers in **Expanded Polystyrene** are designed to maximize temperature control and withstand the hazards of transport. Containers include dome-style lids to maximise shipment capacity, moulded handles, rounded corners and edges, as well as recessed strapping grooves.

Description	Dimensions (hxwxh) (mm)	Qty/Pk	Cat. No.
Shipping Container	152 x 101 x 101	27	216-0202
Shipping Container With Outer Carton	159 x 108 x 140	24	216-0203
Shipping Container	203 x 152 x 101	12	216-0205
Shipping Container With Outer Carton	203 x 152 x 108	12	216-0206
Shipping Container With Outer Carton	203 x 149 x 178	12	216-0207
Shipping Container	203 x 152 x 203	12	216-0208
Shipping Container	203 x 152 x 228	12	216-0209
Shipping Container With Outer Carton	203 x 149 x 255	8	216-0210
Shipping Container	203 x 152 x 304	8	216-0204
Shipping Container	279 x 228 x 228	4	216-0198
Shipping Container With Outer Carton	298 x 254 x 219	6	216-0199
Shipping Container	355 x 355 x 279	4	216-0200
Shipping Container	533 x 335 x 279	2	216-0201



Your European Distribution Partner

Austria

VWR International GmbH
Graumanngasse 7
1150 Wien
Tel.: 01 97 002 0
Fax: 01 97 002 600
E-mail: info@at.vwr.com

Belgium

VWR International bvba/sprl
Haasrode Researchpark Zone 3
Geldenaaksebaan 464
3001 Leuven
Tel.: 016 385 011
Fax: 016 385 385
E-mail:
customerservice@be.vwr.com

Denmark

VWR International ApS
Valhøjs Alle 174-176
2610 Rødovre
Tel.: 043 86 87 88
Fax: 043 86 87 90
E-mail: info@dk.vwr.com

Finland

VWR International Oy
Pihatörmä 1 C 1
02240 Espoo
Tel.: 09 80 45 51
Fax: 09 80 45 52 00
E-mail: info@fi.vwr.com

France

VWR International S.A.S.
Le Périgares – Bâtiment B
201, rue Carnot
94126 Fontenay-sous-Bois cedex
Tel.: 0 825 02 30 30 (0,15 EUR TTC/min)
Fax: 0 825 02 30 35 (0,15 EUR TTC/min)
E-mail: info@fr.vwr.com

Germany

VWR International GmbH
Hilpertstrasse 20a
D - 64295 Darmstadt
Tel.: 0180 570 20 00
Fax: 0180 570 22 22
E-mail: info@de.vwr.com

Ireland

AGB Scientific Ltd - A VWR
International Company
Slaney Close
Dublin Industrial Estate
Dublin 11
Tel: 01 88 22 222
Fax: 01 88 22 333
E-mail info@agb.ie

Italy

VWR International s.r.l.
Via Stephenson 94
20157 Milano (MI)
Tel.: 02 332 03 11
Fax: 800 152 999
E-mail: info@it.vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway

VWR International AS
Kakkelvnskroken 1
P.B. 45, Kalbakken
0901 Oslo
Tel.: 0 2290
Fax: 0 22 90 00 40
E-mail: info@no.vwr.com

Portugal

VWR International – Material de
Laboratorio, Lda. Portugal
Rua Alfredo da Silva 3-C
1300-040 Lisboa
Tel: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Spain

VWR International Eurolab S.L.
Apartado 48
08100 Mollet del Vallés - Barcelona
Tel.: 93 565 54 54
Fax: 93 544 00 03
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: info@se.vwr.com

Switzerland

VWR International AG
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard
Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com